



Hose Solutions

- High Purity TFE Hose Assemblies
- BioPharm USP Class VI Silicones
- USDA Food Grade Rubber
- Industrial Plastics and Rubber
- Validation Documentation and Certifications
- Proudly Made in the U.S.A.



Visit us at www.sentinelprocess.com
or Call 800.345.FLOW Today!

Contact Sentinel Process Systems, Inc. to help with all your Process Component, Instrumentation, and Critical Utility needs.

Hose Solutions



for tough process applications

If your process requires a TFE lined hose, Sentinel can quickly fabricate custom lengths with encapsulated PFA Teflon, Hastelloy or AL6XN Tri-Clamps, Flanges, or Cam-Locks.

Sentinel uniquely serializes every hose assembly we fabricate in order to help you manage your production hose inventory. We also offer a variety of custom tagging and labeling options as well. Our engineers are ready to help you with all your flex hose applications.



PFA Lined
CAM Lock Adaptors



Molded Silicone Labels



Laser Etch Serializing



HIGH PURITY HOSE ASSEMBLIES

- Smooth Bore PTFE/Stainless Steel Braid Silicone outer cover (TSC)
- Rubber Covered FEP Lined
- Stainless Steel or Polypro Braid Covered PTFE
 - Smooth Bore
 - Convoluted
- End Connections:
 - Tri-Clamp
 - ANSI Flanged
 - Cam-Lock
 - MNPT or FNPT

BIOPHARM SILICONE

- Clear
- Nomex Braid Reinforced
- 4-Ply Wire Reinforced
- Convoluted Outer Cover
- Tri-Clamp Ends
 - 15 Ra Finish
 - Electropolished



USDA FOOD GRADE RUBBER

- EPDM
- Buna
- Viton

INDUSTRIAL PLASTICS

- PVC
- PolyPro

SPECIAL OPTIONS

- Color Coding
- Vulcanized Tagging
- Electro Etched
- Custom Labels
- Heat Traced Hose

CUSTOM HOSE RACK

- Polished Stainless Steel
- Heavy Gauge, Rugged Design
- Multiple Mounting Options
- Safely Store Clamps & Gaskets

CORPORATE OFFICE:

3265 Sunset Lane, Hatboro, PA 19040 • P: 215.675.5700 • F: 215.675.9540

NEW ENGLAND:

50 Brigham Street, Marlboro, MA 01752 • P: 508.624.5577 • F: 508.624.7755