

Single-Use Media Storage

Sentinel Process offers a wide range of Media Storage solutions used in single-use applications. Our solutions are designed to meet complex demands for material integrity throughout the process. Available in standard designs or customized solutions to meet specific customer requirements with short lead times.

- **Variety of storage designs to address every application**
- **Various materials including LDPE, EVOH, PTEG, PP, HDPE, PC and Fluoropolymer**
- **Silicone and TPE tubing with numerous termination options**
- **Gamma irradiated between 25 and 40kGy**
- **Assembled in ISO 7 Clean Room**
- **All product contact materials USP Class VI Tested**
- **Full traceability for all materials**
- **Certificates of analysis provided for each lot**
- **Extractables and leachables (E&L) data available for all contact materials**



Visit us at www.sentinelprocess.com or Call 800.345.FLOW Today!

Single-Use Media Storage



2D Bioprocess Bags

50L to 200L
1 to 3 ports, Single, 2 and 3-ply



3D Bioprocess Bags

50L to 1000L
1 to 5 ports, 3-ply



Tank Liners

50L to 1500L
Open or Closed Top, 3-ply



Motion Bioreactor Bags

5L to 50L, 4, 5 and 7 ports,
Single, 2 and 3-ply



Mixing System Bags

50L to 200L
5 ports, 3-ply



Microcarrier Filter Bags

500ml to 4L
Table-Top or Hanging, 2-ply



Cryopreservation Bags

50ml to 20L, 2 or 3 port,
-85°C to 40°C, (-121°F to 104°F)



Rigid Containment Barrels

50L to 200L, Standard, Mixing and
Temperature Controlled



Bag Manifolds

Custom Sizes, Lengths
and Configurations



Media Bottles

60ml to 2L
PTEG



Carboys

2.5L to 40L
Various Materials Available



Also Check
Out Our
Tubing
Solutions
Flyer



CORPORATE OFFICE:

3265 Sunset Lane, Hatboro, PA 19040 • P: 215.675.5700 • F: 215.675.9540

NEW ENGLAND:

50 Brigham Street, Marlboro, MA 01752 • P: 508.624.5577 • F: 508.624.7755