

# Single-Use Tubing Solutions

Sentinel Process offers a broad product portfolio of Tubing solutions used in single-use applications. Our solutions are designed to meet complex demands for material integrity throughout the process. Available in customized designs, sizes and packaging to meet specific customer requirements with short lead times.

- **Tubing set assemblies**
- **Bulk tubing in rolls or spools**
- **Molded manifold assemblies and silicone hygienic ends**
- **Sanitary and Bioprocessing Fittings**
- **Variety of material including high purity silicone and weldable TPE**
- **Gamma irradiated between 25 and 40kGy**
- **Assembled in ISO Class 7 Clean Room**
- **All product contact materials USP Class VI Tested**
- **Full traceability for all materials**
- **Certificates of analysis provided for each lot**
- **Extractables and leachables (E&L) data available for all contact materials**




**Visit us at [www.sentinelprocess.com](http://www.sentinelprocess.com) or Call 800.345.FLOW Today!**

# Single-Use Tubing Solutions




**Tubing Set Assemblies**  
Customized, preassembled  
Silicone and TPE tubing



**Bulk Tubing — Silicone**  
1/32" to 1" ID,  
unreinforced or braided



**Bulk Tubing — TPE**  
1/8" to 1" ID, suitable  
for peristaltic pumps



**Molded Manifolds**  
Silicone and TPE assemblies,  
with molded shapes



**Molded Silicone Ends**  
Smooth bore, eliminates leaks  
in aseptic transfer



**Container Closure Assemblies**  
Stoppers, caps, inserts  
with or without tubing



**Sanitary Fittings**  
Barbed adapters, Reducers, TC's




**Bioprocessing Fittings**  
Quick connects, aseptic connectors,  
Luer locks



**Tube Sealers/Welders**  
Sealers: 0.125" – 0.375" OD  
Welder: 0.375" – 0.75" OD



**Pinch Valves**  
2 body sizes handles  
60 tube sizes



**Ultrasonic Flow Meters**  
Highly accurate, non-invasive  
clamp-on flow sensor



**Also Check  
Out Our  
Media  
Storage  
Solution**



**CORPORATE OFFICE:**

3265 Sunset Lane, Hatboro, PA 19040 • P: 215.675.5700 • F: 215.675.9540

**NEW ENGLAND:**

50 Brigham Street, Marlboro, MA 01752 • P: 508.624.5577 • F: 508.624.7755